


BV Surgical Instrument (356400)



Surgical Kit Instruments

	Ø2.0 Twist Drill	Ø3.0 Twist Drill	F3.5 Taper Drill	F4.0 Taper Drill	F4.5 Taper Drill	F5.0 Taper Drill
L	7mm	351080		351700	351760	351820
	8.5mm	351110		351710	351770	351830
	10mm	351140	351180	351720	351780	351840
	11.5mm	351170	351300	351730	351790	351850
	13mm	351250		351740	351800	351860
	15mm	351380		351750	351810	351870

Lance Drill - Guide Drill	Sidecut Drill	Cortical Drill	Drill Extension	Parallel Pin	Hand Driver
---------------------------	---------------	----------------	-----------------	--------------	-------------

351030 (Ø2.0x13)	351050 (Ø2.0x15)	351610 F3.5 351620 F4.0 351630 F4.5 351640 F5.0	351940 (Ø4.5x26.5)	352010 Ø4.0 352020 Ø5.0	353010 1.2S 353020 1.2L
------------------	------------------	----------------------------------------------------------	--------------------	----------------------------	----------------------------

Torque Driver	Mount Driver	Mount Extension	NoMount Driver	Fixture Driver	Removal Tool
---------------	--------------	-----------------	----------------	----------------	--------------


353030 1.2S 353040 1.2L	354010 S 354020 L	354030 S 354040 L	355050 Mini 355060 Regular	355030 Mini 355040 Regular	355010 Mini 355020 Regular
----------------------------	----------------------	----------------------	-------------------------------	-------------------------------	-------------------------------

Torque Wrench	Depth Gauge	Ratchet Wrench	Service Instrument for Ratchet
---------------	-------------	----------------	--------------------------------


355110	352030	356010	356030
--------	--------	--------	--------

Surgery Procedure

Mini



Regular



Directions for use

- 1) It is an operation guideline for implanting implants in the up and down alveolar bone.
- 2) Using a lance drill pierce the cortical bone and determine the implant position.
- 3) Use the same length twist drill as the implant code. Accompany the drilling with irrigation and pumping in order to keep the heat down from friction.
- 4) After every twist drilling, use depth gauge check the hole depth and floor condition, use parallel pin check orientation.
- 5) Processing exclusive use taper drill for implant diameter and length (Ø3.5 implant: F3.5, Ø4.0 implant: F3.5→F4.0, Ø4.5 implant: F3.5→F4.0→F4.5, Ø5.0 implant: F3.5→F4.5→F5.0, Ø6.0 implant: F6.0, Ø7.0 implant: F6.0→F7.0).
- 6) Use taper cortical drills for implant diameter and length after formation of final drill hole in case of more than hard bone.
- 7) Finally, the fixture is implanted.

Cautions

- When drilling, move the hand piece perpendicularly up and down in a pumping motion.
- To reduce the friction during drilling, provide ample cooling with pre-cooled (5 °C, 41 °F) sterile saline solution.
- The drilling speed must be maintained at 800—1000rpm, the Cortical drilling speed must be maintained at 300—400rpm, and the hole must be created in advance. Dental implant placement should be accomplished at very low speed (15 rpm) or manually.
- Recommended number of use: drill <50 times.
- Recommended number of use: driver <50 times, the maximum allowable torque of 1.2hex driver is 35Ncm.

Dental Implant System

Add: No.26, Huashan Middle Road, Xinbei Zone, Changzhou, Jiangsu, China, 213022 www.bioconcept.cn

Tel:+86-519-67895691 E-mail:intlsales@bioconcept.cn

