

Bioconcept

Dental Implant Systems

Benefit more people



Benefit more people



Bioconcept

Leading Brand for Premium Dental Implants

Bioconcept is a world-class manufacturer of high quality, high-precision dental implants, prosthetics and surgical instruments for the global market. We offer two primary implant systems: **Bioconcept BV** with Tapered Bone Level implants, **Bioconcept BC** with Tissue Level (T-L), Bone Level (B-L), and Bone Level Tapered (B-L-T) implants. Bioconcept offerings are sold in 60+ countries, including USA, Canada, Germany, Italy, France, Australia and China.

We have fully integrated R&D, Production, QA/QC, Sales and Service Support, and utilize only the highest-quality raw materials from leading global suppliers (ZAPP and Sandvik). Our senior management team averages 25+ years of medical device expertise, previously having worked at leading MNCs, such as Medtronic, Johnson & Johnson, Baxter and Komatsu.

At Bioconcept our mission remains: "Benefit more people"



S-L-A

Surface Treatment

- Gold Standard of Surface Treatment
- Promotes bi-directional osseointegration
- Reduces the incidence rate of peri-implantitis

Bioconcept "At-a-Glance"

- Implanted over 1,000,000 across the world
- Our products are sold in 60+ countries
- 25,000 m² state-of-the-art production facilities
- 25+ years for senior management experience
- Premium 'one-stop' solutions since 2007

Key Product Features & Innovation

- Grade 4 (h6) Titanium and Medical Grade Stainless Steel
- Outstanding Surface Treatment with proprietary S.L.A
- Best-in-Class Surgical Kits & Instruments
- Bi-Thread Center Screw
- Strong Self-Retaining Performance
- No Loosening of Screw Structure
- Patented Plug-in Carrier

Multi-Unit Abutments



Equator Prosthetic System



CAD/CAM Abutments



Standard, Mini & Specialty Surgical Kits

Our Premium 'One-Stop' Solutions

Dental Implants

Bioconcept BC

Key BC Features:

Ø (mm): 3.3 | 4.1 | 4.8

L (mm): 6 | 8 | 10 | 12 | 14 | 16


Ti Grade 4 (h6) (ZAPP) and S-L-A Surface

Tissue Level (T-L)

- Connection: Octagon 16° morse-taper
- Type: S, SP, TE (NN, RN, WN)

Bone Level (B-L) and Tapered (B-L-T)

- Connection: Cross (NC, RC)
- B-L-T: Apically at 0.8 mm pitch



Bioconcept BV

Key BV Features:

Ø (mm): 3.5 | 4.0 | 4.5 | 5.0 | 6.0 | 7.0

L (mm): 7 | 8.5 | 10 | 11.5 | 13 | 15

Ti Grade 4 (h6) (ZAPP) and S-L-A Surface

Tapered Bone Level

- Connection: Hexagon, 11° morse-taper
- Platform switching and open thread
- Double-lead thread
- Bi-thread center screw innovation
- Plug-in patented carrier



Prosthetics

Full-Line of Prosthetic Solutions

Single Tooth
Multi-Units
Full Denture
Full Overdenture

BC: T-L



Solid | Minimal Height |
Octagonal (Straight & Angled) | Equator

BC: B-L/B-L-T



Anatomic | Meso | Multi-Function
(Straight & Angled) | Equator

BV:



Transfer | Angled | Multi-Function
(Straight & Angled) | Equator

Surgical Kits

Best-in-Class Surgical Kits & Instruments



Surgical Kits



BV Guided



Mini Kits



Specialty Kits



- Implants
- Prosthetic
- Rescue
- Bone Expansion
- Bone Collector
- Sinus
- Osteo
- Surgery

Digital Solutions

Tibase



- Extended Prosthetics Solution
- Maximum Design Freedom

Indications: Single | Multi-

Fixation: Cement | Screw

Workflow: Traditional | Digital

Implant System: BC | BV

Platforms: NC | RC | RN | WN

Spiral groove: Straight

CAD/CAM Ti-Customized



Pre-Milled Abutments

- BC
- BV
- Dentium
- Nobel
- ICX
- MG
- Astra

- Highest Accuracy Guaranteed
- Best Solution for Aesthetics and LT Stability

Benefit more people

Bioconcept

State-of-the-Art Production & QA/QC

In 2023, Bioconcept opened its state-of-the-art 25,000 m² production facility in Changzhou, China, establishing one of the largest and most advanced facilities in the world.

Annual Capacity: 800,000+ Implants
(Reaching 3.0 million by 2026)


- 100+ Swiss machines operated continuously (reaching 250 by 2025)
- 1,500 m² clean workshop Class 8 and ISO 13485 certified
- Premier X-ray sterilization with Cobalt 60 irradiation systems
- World-class QA/QC (3D imaging, dynamic photocast, SEM, etc.)



Our Global Network

Bioconcept is recognized by customers across the globe for our highest-quality guaranteed products, our professional and timely service, and our educational and clinical support network.

Our Products Are Sold In
More Than 60 Countries



100+ sales and service team members



350+ courses; 4k+ dentists trained per year



70+ online trainings per year with 30k followers



30+ academic conferences per year



40+ global dental exhibitions per year

Contact Us

Web: www.bioconcept.com

Email: support@bioconcept.com

Tel: +86-519-67895691

Add: China Headquarters: No.1019, Liaohe Road, Xinbei District, 213000, Changzhou, Jiangsu, PRC

North America Office: Suite 405 - 1250 Homer Street, Vancouver, BC, V6B 2Y5, Canada



Web