




Mini



Regular



		Mini				Regular								
S-L-A Implant	L	Mount Φ3.5	NoMount Φ3.5	Mount Φ4.0	NoMount Φ4.0	Mount Φ4.5	NoMount Φ4.5	Mount Φ5.0	NoMount Φ5.0	Mount Φ6.0	NoMount Φ6.0	Mount Φ7.0	NoMount Φ7.0	
		7mm	311010	311010G	312010	312010G	313010	313010G	314020	314020G	315020	315020G	316020	316020G
8.5mm	311020	311020G	312020	312020G	313020	313020G	314030	314030G	315030	315030G	316030	316030G		
10mm	311030	311030G	312030	312030G	313030	313030G	314040	314040G	315040	315040G	316040	316040G		
11.5mm	311040	311040G	312040	312040G	313040	313040G	314050	314050G	315050	315050G	316050	316050G		
13mm	311050	311050G	312050	312050G	313050	313050G	314060	314060G	315060	315060G	316060	316060G		
15mm			312060	312060G	313060	313060G	314070	314070G						
Lance Drill														
Pilot Drill	L													
		7mm	351080	351080	351080	351080	351080	351080	351080	7mm	352510	352510	352510	
		8.5mm	351110	351110	351110	351110	351110	351110	351110	8.5mm	352510	352510	352510	
		10mm	351140	351140	351140	351140	351140	351140	351140	10mm	352510	352510	352510	
		11.5mm	351170	351170	351170	351170	351170	351170	351170	11.5mm	352520	352520	352520	
		13mm	351250	351250	351250	351250	351250	351250	351250	13mm	352520	352520	352520	
Twist Drill	L													
		7mm	351180/351300	351700	351700	351700	351760	351700	351820	352530	352530	352590	352590	
		8.5mm		351710	351710	351770	351760	351710	351820	352540	352540	352600	352600	
		10mm		351720	351720	351780	351780	351720	351840	352550	352550	352610	352610	
		11.5mm		351730	351730	351790	351790	351730	351850	352560	352560	352620	352620	
		13mm		351740	351740	351800	351800	351740	351860	352570	352570	352630	352630	
	L	15mm	351300	351750	351750	351810	351750	351870	351750	351870				
		7mm												
		8.5mm												
		10mm												
		11.5mm												
		13mm												
Taper Cortical Drill														
1.Cover Screw														
2.Healing Abutment		323010(H3.0) 323020(H4.0) 323030(H5.0) 323040(H7.0)				322010(H0.4) 324010(H3.0) 324020(H4.0) 324030(H5.0) 324040(H7.0)				324510(H3.0) 324520(H4.0) 324530(H5.0) 324540(H7.0)				
		323110(H3.0) 323120(H4.0) 323130(H5.0) 323140(H7.0)				324110(H3.0) 324120(H4.0) 324130(H5.0) 324140(H7.0)				324210(H3.0) 324220(H4.0) 324230(H5.0) 324240(H7.0)				
		323310(H3.0) 323320(H5.0) 323330(H7.0)				324310(H3.0) 324320(H4.0) 324330(H5.0) 324340(H7.0)				324410(H3.0) 324420(H4.0) 324430(H5.0) 324440(H7.0)				

Bioconcept

BV

Tapered Bone Level Implant System


Tel: +86-519-67895691
 E-mail: intlsales@bioconcept.cn
 https://www.bioconcept.cn
 Add: No.1019, Liaohe Road, Xinbei District, Changzhou, Jiangsu, 213125, China

Temporary Abutment /Esthetic-low Healing Cap

Temporary Abutment

H	GH	Type	D	Φ4.0
10	1.0	Hex	371010	
		Non-hex	371020	
	3.0	Hex	371030	
		Non-hex	371040	


* Tightening torque: 20-30Ncm



Temporary Abutment

H	GH	Type	D	Φ4.5
10	1.0	Hex	371050	
		Non-hex	371060	
	3.0	Hex	371070	
		Non-hex	371080	

* Tightening torque: 20-35Ncm



Healing Cap

Equator abutment number	325010	325020	325030
	338010	338050	338250
	338020	338060	338260
	338030	338230	
	338040	338240	
	338210		
	338220		

Impression Component

Impression Coping

Manual Impression

L	D	Φ4.0
11		361110
15		361120
Open tray		
		361160
		361170
Closed tray		

Digital Impression

L	D	Φ3.8
10		361310
15		361320
Scanbody		

Manual Impression

Φ4.8
362470
Open tray
Φ4.8
362480
Closed tray

Digital Impression

Φ4.8
362640
Scanbody

Manual Impression

L	D	Φ4.0	Φ5.0	Φ6.0
11		362210	362220	362230
15		362240	362250	362260
Open tray				
		362300	362320	362340
		362310	362330	362350
Closed tray				

Digital Impression

L	D	Φ4.5
10		362570
15		362580
Scanbody		

Manual Impression

Φ4.8
362470
Open tray
Φ4.8
362480
Closed tray

Digital Impression

Φ4.8
362640
Scanbody

Analog



Prosthetic Restoration

Cement Retained Restoration

Cement Retained Restoration

Screw Retained Restoration

Cement Retained Restoration

Cement Retained Restoration

Screw Retained Restoration

Retentive Anchor

Transfer Abutment

H	GH	D	Φ4.5
5.5	1.0	Hex	331010
		Non-hex	331020
	3.0	Hex	331030
		Non-hex	331040
		Hex	331050
7.0	1.0	Hex	331510
		Non-hex	331520
	3.0	Hex	331530
		Non-hex	331540
		Hex	331550

* Tightening torque: 30Ncm

Angled Abutment

H	GH	Type	D	Φ4.5
8.0	2.0	A	333010	
		B	333020	
	4.0	A	333040	
		B	333050	

* Tightening torque: 30Ncm

Multi Abutment

GH	D	Φ4.8
1.0		337010
2.0		337020
3.0		337030
4.0		337040
5.0		337050

* Tightening torque: 30Ncm

Multi-Angled Abutment

A	GH	D	Φ4.8
17°	2.5		337130
	3.0		337140
	4.0		337150
30°	3.5		337160
	4.0		337170
	5.0		337180

* Tightening torque: 25Ncm

Transfer Abutment

H	GH	D	Hex	Non-Hex
5.5	1.0	Φ4.5	331060	331160
		Φ5.0	331170	331320
	3.0	Hex	331080	331180
		Non-hex	331190	331340
		Hex	331300	331440
7.0	1.0	Φ4.5	331110	331210
		Φ5.0	331220	331360
	3.0	Hex	331120	331220
		Non-hex	331230	331380
		Hex	331130	331240
* Tightening torque: 35Ncm				

Angled Abutment

H	GH	Type	Φ5.0	Φ6.0
8.0	2.0	A	333070	333130
		B	333080	333140
	4.0	A	333100	333160
		B	333110	333170

* Tightening torque: 35Ncm

Multi Abutment

GH	D	Φ4.8
1.0		337070
2.0		337080
3.0		337090
4.0		337100
5.0		337110

* Tightening torque: 30Ncm

Multi-Angled Abutment

A	GH	D	Φ4.8
17°	2.5		337190
	3.0		337200
	4.0		337210
	3.5		337220
30°	4.0		337230
	4.0		337240
	5.0		337250

* Tightening torque: 30Ncm

Equator Abutment

Gingiva Height	Mini	Regular
1	338010	338210
2	338020	338220
3	338030	338230
4	338040	338240
5	338050	338250
6	338060	338260

* Tightening torque: 35Ncm

Titanium Coping

D	H	GH	Hex	Non-Hex
Φ4.0	3.0	1.0	813010	813010N
		2.0	813020	813020N
	5.0	1.0	813030	813030N
		2.0	813040	813040N

* Tightening torque: 30Ncm

Esthetic-low Temporary Cylinder

377010

* Tightening torque: 20Ncm

D	H	GH	Hex	Non-Hex
Φ4.5	3.0	1.0	813050	813050N
		2.0	813060	813060N
	5.0	1.0	813070	813070N
		2.0	813080	813080N

* Tightening torque: 35Ncm

Esthetic-low Temporary Cylinder

377010

* Tightening torque: 20Ncm

SmartBox Male Processing Package 378100

Equator Processing Male	Equator Block-out Spacers	Equator Replacement Male	
		2.7kg	0.6kg
		378020v4	378050v4
		378030v4	378040v4

TI Screw



Auxiliary Parts Screw

Instrument

Hex Hand Drill

353010 (Short)
353020 (Long)

Hex Torque Drill

353030 (Short)
353040 (Long)

Multi Abutment Outer Driver

354050

Multi-angled Abutment Carrier

353060

Hex Hand Drill

353010 (Short)
353020 (Long)

Hex Torque Drill

353030 (Short)
353040 (Long)

Multi Abutment Outer Driver

354050

Multi-angled Abutment Carrier

353050

BV Equator Square Screwdriver for ratchet

358010

Esthetic-low Polishing Protector

357020

Esthetic-low Polishing Protector

357020

Insertion-Extraction Tool

358020