







Bioconcept BC Tissue Level Implant System Tel: +86-519-67895691 E-mail: intlsales@bioconcept.cn https://www.bioconcept.cn Add: No.1019, Liaohe Road, Xinbei District, Changzhou, Jiangsu, 213125, China		Φ3.3mm				Φ4.1mm			Φ4.8mm					
		Standard Plus NN Φ3.5mm  Φ3.3mm	Standard RN Φ4.8mm  Φ3.3mm	Standard Plus RN Φ4.8mm  Φ3.3mm	Tapered RN Φ4.8mm  Φ3.3mm	Standard RN Φ4.8mm  Φ4.1mm	Standard Plus RN Φ4.8mm  Φ4.1mm	Tapered RN Φ4.8mm  Φ4.1mm	Standard RN Φ4.8mm  Φ4.8mm	Standard Plus RN Φ4.8mm  Φ4.8mm	Standard WN Φ6.5mm  Φ4.8mm	Standard Plus WN Φ6.5mm  Φ4.8mm	Tapered WN Φ6.5mm  Φ4.8mm	
		SPΦ3.3NN	SΦ3.3RN	SPΦ3.3RN	TEΦ3.3RN	SΦ4.1RN	SPΦ4.1RN	TEΦ4.1RN	SΦ4.8RN	SPΦ4.8RN	SΦ4.8WN	SPΦ4.8WN	TEΦ4.8WN	
S-L-A Implants	6 mm					013060	012050		013120	012100	013170	012150		6 mm
	8 mm	011010	013010	012010	014010	013070	012060	014050	013130	012110	013180	012160		8 mm
	10mm	011020	013020	012020	014020	013080	012070	014060	013140	012120	013190	012170	014090	10mm
	12mm	011030	013030	012030	014030	013090	012080	014070	013150	012130	013200	012180	014100	12mm
	14mm	011040	013040	012040	014040	013100	012090	014080	013160	012140			014110	14mm
16mm			013050			013110								16mm
Round Burr		 Round Burr Φ1.4 051010  Round Burr Φ2.3 051020  Round Burr Φ3.1 051030				 Round Burr Φ1.4 051010  Round Burr Φ2.3 051020  Round Burr Φ3.1 051030			 Round Burr Φ1.4 051010  Round Burr Φ2.3 051020  Round Burr Φ3.1 051030					
Drill		 Pilot Drill Φ2.2 051040(short) 051050(long)  Twist Drill-1 Φ2.8 051060(short) 051070(long)				 Pilot Drill Φ2.2 051040(short) 051050(long)  Twist Drill-1 Φ2.8 051060(short) 051070(long)  Twist Drill-2 Φ3.5 051080(short) 051090(long)			 Pilot Drill Φ2.2 051040(short) 051050(long)  Twist Drill-1 Φ2.8 051060(short) 051070(long)  Twist Drill-2 Φ3.5 051080(short) 051090(long)  Twist Drill-3 Φ4.2 051100(short) 051110(long)					
Depth Gauge		 Depth Gauge (Alignment Pin) Φ2.2 053010  Depth Gauge Φ2.2/Φ2.8 053020				 Depth Gauge (Alignment Pin) Φ2.2 053010  Depth Gauge Φ2.2/Φ2.8 053020  Depth Gauge Φ3.5 053030			 Depth Gauge (Alignment Pin) Φ2.2 053010  Depth Gauge Φ2.2/Φ2.8 053020  Depth Gauge Φ3.5 053030  Depth Gauge Φ4.2 053040					
Profile Drill		 SP Profile Drill 051120(short) 051130(long)  TE Profile Drill 051180(short) 051190(long)				 SP Profile Drill 051140(short) 051150(long)  TE Profile Drill 051200(short) 051210(long)			 SP Profile Drill 051160(short) 051170(long)  TE Profile Drill 051220(short) 051230(long)					
Tap		 Tap Φ3.3 051240(short)  Tap Φ3.3 051250(long)  Adapter Tap Φ3.3 051260				 Tap Φ4.1 051270(short)  Tap Φ4.1 051280(long)  Adapter Tap Φ4.1 051290			 Tap Φ4.8 051300(short)  Tap Φ4.8 051310(long)  Adapter Tap Φ4.8 051320					
Healing Cap		NN  H3.4mm 021020  H3.4mm 021030		RN  Closure Screw 022010  H1.5mm 022020  H2.0mm 022030  H3.0mm 022040  H4.5mm 022050  H2.0mm 022060  H3.5mm 022070						WN  Closuer Screw 023010  H2.0mm 023020  H3.0mm 023030  H4.5mm 023040  H2.0mm 023050				
		NN=Narrow Neck		RN=Regular Neck		WN=Wide Neck		SP=Standard Plus	SP=Standard	TE=Tapered Effect		S-L-A=Sand-blasted, Large grit, Acid-etched		KL31

Bioconcept BC Tissue Level Implant System Tel: +86-519-67895691 E-mail: intlsales@bioconcept.cn https://www.bioconcept.cn Add: No.1019, Liaohe Road, Xinbei District, Changzhou, Jiangsu, 213125, China	NN(Narrow Neck)	RN(Regular Neck)						WN(Wide Neck)				Overdenture																					
	Φ3.5mm 	Φ4.8mm 						Φ6.5mm 				Abutment Level Impression																					
	Octagonal System	Octagonal System				Solid System	Hybrid	Octagonal System		Solid System	Equator																						
Transfer Parts	NN Impression Cap with integral guide screw 061010 NN "Snap-on" Impression cap 061020 	Manual Impression RN Impression Cap 062170 RN Impression Cap 062180  Digital Impression RN Scanbody 860010 Digital Analog 062220 RN Positioning Cylinder, red 062050 RN Impression Cap built-in handle 062030 RN Impression Cap screwed 062040 RN Impression Cap 062010 RN Impression Cap 062160 				RN Positioning Cylinder, yellow 062060 RN Positioning Cylinder, grey 062070 RN Positioning Cylinder, blue 062080  Analog for RN Solid Abutment H4mm, yellow 062100 H5.5mm, grey 062110 H7mm, blue 062120		Manual Impression WN Impression Cap 063130 WN Impression Cap 063140  Digital Impression WN Scanbody 860020 Digital Analog 063180 WN Positioning Cylinder, grey 063030 WN Impression Cap with integral guide screw 063020 WN Impression Cap 063010 WN Impression Cap 063120 	WN Positioning Cylinder, green 063040 WN Positioning Cylinder, brown 063050  Analog for WN Solid Abutment H4mm, green 063070 H4mm, brown 063080	368010 368020 Retentive Anchor Equator Abutment  <table><tr><td></td><td>GH</td><td>RN</td><td>WN</td></tr><tr><td rowspan="6">Gingiva Height</td><td>1</td><td>034010</td><td>034210</td></tr><tr><td>2</td><td>034020</td><td>034220</td></tr><tr><td>3</td><td>034030</td><td>034230</td></tr><tr><td>4</td><td>034040</td><td>034240</td></tr><tr><td>5</td><td>034050</td><td>034250</td></tr><tr><td>6</td><td>034060</td><td>034260</td></tr></table> * Tightening torque: 35Ncm		GH	RN	WN	Gingiva Height	1	034010	034210	2	034020	034220	3	034030	034230	4	034040	034240	5	034050	034250	6	034060	034260
		GH	RN	WN																													
Gingiva Height	1	034010	034210																														
	2	034020	034220																														
	3	034030	034230																														
	4	034040	034240																														
	5	034050	034250																														
	6	034060	034260																														
Temporary Prosthetics		H9mm for crown 072010 H9mm for bridge 072010N * Tightening torque:15-35Ncm					H9mm for crown 073010 H9mm for bridge 073010N * Tightening torque:15-35Ncm																										
Final Prosthetics	Cement-retained	Screw-retained	Screw-retained or cement-retained		Cement-retained	Cement-retained	Bar	Retentive Anchor	Screw-retained	Cement-retained	Cement-retained	Cement-retained																					
Abutments Titanium			     			  	 			 		 																					
TiBase Abutment	can be trimmed 031010 15° Angled 031020 20° Angled 031030	H4mm Crown 810010 H6mm Crown 810020 H3.5mm Bridge/bar 810010N H4mm Bridge/bar Cylindrical 811150N							TiBase Abutment H4.5mm Crown 810030 H6.5mm Crown 810040 H4.8mm Bridge/bar 810030N H3.5mm Bridge/bar Cylindrical 811160N																								
Copings Plastic		Bridge/Bar 090030 Crown 090040	Shoulder For RN Angled 090080 For Angled RN Abutment 090070		Bridge 090050 Crown 090060	Bridge 090010 Crown 090020			Bridge 090110 Crown 090120	Shoulder For WN15° Angled Abutments 090150	Bridge 090130 Crown 090140	Bridge 090090 Crown 090100																					
Auxiliary Parts Screws	 * Tightening torque:35Ncm	 * Tightening torque:15Ncm							 * Tightening torque:15Ncm			Equator Insertion-Extraction Tool  358020																					