

Equator Abutment

Overdenture solution

- Smallest Attachment on the market
- Better Aesthetic Result
- Severe Implant Divergency Solution



Smart box

- Correct Divergency up to 50°



Low profile titanium implant abutment

- Gingival Heights Available: 1 mm up to 6 mm
- Available on all the implant platforms and connections

Retentive caps

- Different elastic retention degrees
- Different diameters
- Allows for correction of implant divergency
- Shock absorbing system
- Trauma free insertion of prosthesis














Equator Abutment



Equator Abutment						
System	BV		BC-TL		BC-BL	
	Mini	Regular	RN	WN	NC	RC
GH						
1mm	338010	338210	034010	034210	—	134210
2mm	338020	338220	034020	034220	134020	134220
3mm	338030	338230	034030	034230	134030	134230
4mm	338040	338240	034040	034240	134040	134240
5mm	338050	338250	034050	034250	134050	134250
6mm	338060	338260	034060	034260	134060	134260

Tightening Torque: 35N·cm

BV Healing caps			
Art.No			
	325010	325020	325030
Equator Abutment No	338010	338050	338250
	338020	338060	338260
	338030	338230	
	338040	338240	
	338210		
	338220		

Component & Tools						
Impression Coping  368010 Analog  368020	SmartBox Male Processing Package 378100 Equator Block-out Spacers Equator Processing Male Equator Block-out Spacers 2.7kg    378020			Equator Replacement Male 1.8kg 1.2kg 0.6kg    378030 378040 378050		
	Screw Driver BC BV   055130 358010	Insertion-Extraction Tool INSERTOR Side  Central holder  EXTRACTOR Side  358020				